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**DESCRIPTION**

OSSIX® PLUS and OSSIX® VOLUMAX are biodegradable and biocompatible collagen membranes intended for use during the process of guided tissue and bone regeneration. The membranes are produced by standardized and controlled manufacturing procedures. The collagen is extracted from porcine tendons subjected to veterinarian inspection and is purified to prevent hypersensitivity reactions. The membrane is packaged in a sealed double blister contained in a cardboard box and is terminally sterilized by ethylene oxide (EtO). The membrane is intended for a single use.

**PROPERTIES**

The membrane has been demonstrated to be biocompatible. Animal and human clinical testing show a low potential to induce hypersensitivity.

The membrane has a porous structure; the size of the pores is small enough to occlude gingival cells but large enough to allow the passage of fluids, nutrients and plasma proteins, which are necessary to support healing (see reference h below).

The membrane is not self-supporting and therefore it is recommended for use in conjunction with support such as an autogenous bone graft, allograft, xenograft, osteoconductive and/or inductive bone substitute, or a mixture of these.

The membrane does not dissolve or disintegrate when wet.

The membrane conforms easily to the shape of the alveolar ridge.

Animal studies shown that the membrane degradation is completed within approximately 6 months.

**INDICATIONS**

The membrane is a resorbable collagen membrane intended for use during the process of guided bone regeneration (GBR) and guided tissue regeneration (GTR) as a biodegradable barrier for:

1. Ridge augmentation for later implant insertions.
2. Simultaneous ridge augmentation and implant insertions.
3. Ridge augmentation around implants inserted in delayed extraction sites.
4. Ridge augmentation around implants inserted in immediate extraction sites.
5. Alveolar ridge preservation consequent to tooth (teeth) extraction(s).
6. Over the window in lateral window sinus elevation procedures.
7. In implants with vertical bone loss due to infection, only in cases where satisfactory debridement and implant surface disinfection can be achieved.
8. In intra bony defects around teeth.
9. For treatment of recession defects, together with coronally positioned flap.
10. In furcation defects in multi rooted teeth.

**CONTRAINDICATIONS**

The membrane must not be used in:

1. Patients with known collagen hypersensitivity.
2. Patients with sensitivity to porcine-derived materials.
3. Patients suffering from autoimmune diseases and connective tissue diseases, such as: systemic lupus erythematosus, dermatomyositis etc.

**WARNINGS AND PRECAUTIONS**

1. The membrane is intended for a single use. Do not reuse or re-sterilize the membrane.
2. Treatment of high risk patients, such as: smokers, patients suffering from uncontrolled diabetes mellitus, and uncontrolled periodontal disease may be impaired.
3. The safety of treatment with the membrane in pregnant and nursing women and in children has not been yet established.
4. The outcome of regenerative procedures may be impaired in patients suffering from untreated periodontitis. Infection control and good oral hygiene should be achieved prior to surgical intervention.

**ADVERSE EVENTS**

1. Clinical and post marketing experience with the membrane reveals an excellent safety profile.
2. Adverse reactions with the membrane were not observed.
3. Yet, as the membrane is of a collagen origin, allergic reactions (e.g. erythema, swelling, induration and/or pruritus at treatment site) may not be entirely excluded.

**DIRECTIONS FOR USE**

1. Special instructions for use in periodontology:  
A basic requirement for successful periodontal treatment includes eradicating the underlying bacterial infection as well as adequate oral hygiene. Therefore, prior to surgical intervention, patients must receive a hygiene phase of treatment, consisting of oral hygiene instructions, scaling and root planing, and occlusal adjustment when indicated. A postoperative maintenance phase can help to ensure long-term therapeutic success.
2. The bony defect should be exposed by full thickness mucoperiosteal flaps.
3. All soft tissues should be removed.
4. In GTR the root surface should be carefully debrided and planed. Root conditioning should be considered.
5. The cortical plate can be perforated in order to allow osteogenic tissues from the bone marrow to colonize the regenerating site.
6. By using sterile a-traumatic instruments and sterile gloves rinsed with sterile saline, the membrane is removed aseptically from the package.
7. The membrane should be immersed (the inner blister can be used as a dish) for 30 seconds in sterile saline, to allow for its expansion to its final dimensions (OSSIX® PLUS: 15x25 mm, 25x30 mm, 30x40 mm; OSSIX® VOLUMAX: 10x12.5 mm, 15x25 mm, 25x30 mm and 10x40 mm). Initial trimming to the estimated final size may be performed prior to sterile saline immersion.
8. Trimming to the required dimensions: it is recommended that the membrane extends 3-4 mm beyond the margins of the defect. One-two mm of uncovered bone to adjacent teeth must be allowed.
9. The membrane is cut with sterile scissors over a sterile container; try-in over the defect with a-traumatic instruments should be performed.
10. The site to be augmented should be filled with a space-maintaining material. The user should follow the manufacturer's instructions for the material used.

11. The membrane should be secured under the lingual flap, then a bone graft placed, and the membrane carefully adapted over the defect. The membrane will adhere to the underlying tissue; additional fixation of the membrane should be considered. Fixation with screws or tacks or suturing through the membrane is not recommended as it may tear the membrane. Fixation with overlying sutures is advised. This may be achieved by anchoring a mattress suture in the apical periosteum buccally and lingually.
12. The mucoperiosteal flaps are sutured while ensuring that the tissue is not under tension. Do not compromise blood supply to the defect area.
13. In GTR, the use of a periodontal dressing may be considered.

**GUIDELINES FOR THE PATIENT**

The success of any surgical treatment depends on fulfilling the directions for use along with guiding the patient, as follows:

1. Preoperative patient's education regarding adequate oral hygiene and meticulous prophylaxis.
2. Postoperative patient's care, e.g.:
  - a. Soft diet, avoidance of contact with tongue, hard food or denture.
  - b. Avoidance of contact with hot temperature food or liquids that may cause early disintegration of the collagen matrix.
  - c. Following suture removal, rinsing with chlorhexidine for one minute twice a day or according to the manufacturer's instructions.

**POSTOPERATIVE REMINDERS**

1. Clinical experience with OSSIX® PLUS reveals no inflammatory signs following accidental exposure. The membrane degrades slowly in the oral environment and the exposed area is covered by connective tissue and epithelium within a few weeks (see references a-h below).
2. Possible complications with any surgery in the oral and maxillofacial region include: infection, flap slough, perforation, abscess formation, bone loss, pain, soft tissue irregularities, and complications associated with the use of anesthesia.
3. Depending upon the type and severity of the complications, as judged by the dental surgeon, membrane removal may be indicated.

**STORAGE AND HANDLING**

1. The membrane should be used by skilled, experienced and/or trained dental surgeons.
2. The material should be handled using sterile gloves or sterile a-traumatic instruments.
3. Placement of the membrane should be performed after membrane's immersion in saline for 30 seconds.
4. Do not use the membrane if it is torn and/or damaged.
5. Do not use the membrane, in the event that the sterile packaging is opened and/or damaged.
6. Any remaining / unused membrane should be discarded according to local regulations.
7. The membrane should be stored at temperatures between 15-30°C (59-86°F).
8. Do not use the membrane after the expiration date.
9. The membrane is intended for a single use. Do not reuse.

**HOW SUPPLIED**

1. The membrane is supplied in a double blister pack, for single use only. Each pack contains one membrane.
2. OSSIX® PLUS is available in three sizes: 15x25 mm, 25x30 mm, 30x40 mm. OSSIX® VOLUMAX is available in four sizes: 10x12.5 mm, 15x25 mm, 25x30 mm and 10x40 mm.

For any further assistance, support or questions, please contact the local distributor or manufacturer.

**REFERENCES**

- a. Moses *et al.* 2005 – Moses O, Pitaru S, Artzi Z, Nemcovsky CE. (2005) Healing of dehiscence-type defects in implants placed together with different barrier membranes: a comparative clinical study. *Clin Oral Implants Res.* 16, 210-219.
- b. Friedmann *et al.* 2001 – Friedmann A, Strietzel FP, Marezki B, Pitaru S, Bernimoulin JP. (2001) Observations on a new collagen barrier membrane in 16 consecutively treated patients. Clinical and histological observations. *J Periodontol.* 72(11):1616-1623.
- c. Friedmann *et al.* 2002 – Friedmann A, Strietzel FP, Marezki B, Pitaru S, Bernimoulin JP. (2002) Histological assessment of augmented jaw bone utilizing a new collagen barrier membrane compared to a standard barrier membrane to protect a granular bone substitute Material. *Clin. Oral Impl. Res.* 13, 587-594.
- d. Klinger *et al.* 2010 – Klinger A, Asad R, Shapira L, Zubery Y. (2010) In vivo degradation of collagen barrier membranes exposed to the oral cavity. *Clin. Oral Impl. Res.* 21, 873-876.
- e. Hoang and Mealey 2012 – Hoang TN, Mealey BL. (2012) Histologic comparison of healing after ridge preservation using human demineralized bone matrix putty with one versus two different – sized bone particles. *J Periodontol.* 83:174-181.
- f. Beilittum *et al.* 2010 – Beilittum I, Artzi Z, Nemcovsky CE. (2010) Clinical evaluation of particulate allogeneic with and without autogenous bone grafts and resorbable collagen membranes for bone augmentation of atrophic alveolar ridges. *Clin. Oral Impl. Res.* 21: 1242-1250.
- g. Scheyer and McGuire 2014 – Scheyer T. and McGuire MK. (2014) Evaluation of Premature Membrane Exposure and Early Healing in Guided Bone Regeneration of Peri-Implant Dehiscence and Fenestration Defects With a Slowly Resorbing Porcine Collagen Ribose Cross-Linked Membrane: A Consecutive Case Series. *Clinical Advances in Periodontics.* 5: 165-170.
- h. Friedmann *et al.* 2008 – Friedmann A, Dehnhardt J, Kleber BM, and Bernimoulin, JP. (2008) Cytocompatibility of collagen and ePTFE membranes on osteoblast-like cells in vitro. *J Biomed Materials Res A* 86:935-941.

**Manufacturer 제조원**

Datum Dental Ltd.  
1 Bat Sheva St.  
P.O.Box 6170  
Lod 7116003  
ISRAEL  
Tel: 972-8-6705400  
Fax: 972-8-6705429  
www.ossidental.com