## OSSIX® PLUS

### Resorptivna kolagena membrana

#### UPUTSTVO ZA UPOTREBU OSSIX® PLUS

##### OPIS

**OSSIX® PLUS** je biorazgradiva i biokompatibilna kolagena membrana namenjena za upotrebu tokom vođene regeneracije tkiva i kostiju. Proizvedena je prema standardizovanim i kontrolisanim procedurama proizvođača.

Kolagen se ekstrahuje iz svinjskih tetiva podvrgnutih veterinarskom pregledu i prečišćava se kako bi sprečile reakcije vezane za preosetljivost.

**OSSIX® PLUS** je pakovan u zapečaćeno dvostruko blister pakovanje, koje se nalazi u kartonskoj kutiji i na kraju sterilisano etilen oksidom (EtO).

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
### REFERENCE

- a. Moses *et al.* 2005 – Moses O, Pitaru S, Artzi Z, Nemcovsky CE. (2005) Healing of dehiscence-type defects in implants placed together with different barrier membranes: a comparative clinical study. *Clin Oral Implants Res.* 16, 210-219.
- b. Friedmann *et al.* 2001 – Friedmann A, Strietzel FP, Marezki B, Pitaru S, Bernimoulin JP. (2001) Observations on a new collagen barrier membrane in 16 consecutively treated patients. Clinical and histological observations. *J Periodontol.* 72(11):1616-1623.
- c. Friedmann *et al.* 2002 – Friedmann A, Strietzel FP, Marezki B, Pitaru S, Bernimoulin JP. (2002) Histological assessment of augmented jaw bone utilizing a new collagen barrier membrane compared to a standard barrier membrane to protect a granular bone substitute Material. *Clin. Oral Impl. Res.* 13, 587–594.
- d. Klinger *et al.* 2010 – Klinger A, Asad R, Shapira L, Zubery Y. (2010) In vivo degradation of collagen barrier membranes exposed to the oral cavity. *Clin. Oral Impl. Res.* 21, 873–876.
- e. Hoang and Mealey 2012 – Hoang TN, Mealey BL. (2012) Histologic comparison of healing after ridge preservation using human demineralized bone matrix putty with one versus two different – sized bone particles. *J Periodontol.* 83:174-181.

- f. Beitlitum *et al.* 2010 – Beitlitum I, Artzi Z, Nemcovsky CE. (2010) Clinical evaluation of particulate allogeneic with and without autogenous bone grafts and resorbable collagen membranes for bone augmentation of atrophic alveolar ridges. Clin. Oral Impl. Res. 21: 1242–1250.
- g. Scheyer and McGuire 2014 – Scheyer T. and. McGuire MK. (2014) Evaluation of Premature Membrane Exposure and Early Healing in Guided Bone Regeneration of Peri-Implant Dehiscence and Fenestration Defects With a Slowly Resorbing Porcine Collagen Ribose Cross-Linked Membrane: A Consecutive Case Series. Clinical Advances in Periodontics. 5: 165-170.
- h. Friedmann *et al.* 2008 – Friedmann A, Dehnhardt J, Kleber BM, and Bernimoulin, JP. (2008) Cytobiocompatibility of collagen and ePTFE membranes on osteoblast-like cells in vitro. J Biomed Materials Res A 86:935-941.

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Datum Dental Ltd.
MedNet EC-REP III GmbH