**OSSIX® VOLUMAX**  
**Resorptivna kolagena membrana**  
**UPUTSTVO ZA UPOTREBU OSSIX® VOLUMAX**

**OPIS**

**OSSIX® VOLUMAX** je biorazgradiva i biokompatibilna kolagena membrana namenjena za upotrebu tokom vođene tkivne i koštane regeneracije. Proizvedena je prema standardizovanim i kontrolisanim procedurama proizvođača.

Kolagen se ekstrahuje iz svinjskih tetiva podvrgnutih veterinarskom pregledu i prečišćava se kako bi sprečile reakcije preosetljivosti.

**OSSIX® VOLUMAX** je pakovan u zapečaćeni dvostruki blister, koji se nalazi u kartonskoj kutiji i na kraju sterilisan etilen oksidom (EtO).

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
### REFERENCE

- a. Moses *et al.* 2005 – Moses O, Pitaru S, Artzi Z, Nemcovsky CE. (2005) Healing of dehiscence-type defects in implants placed together with different barrier membranes: a comparative clinical study. *Clin Oral Implants Res.* 16, 210-219.
- b. Friedmann *et al.* 2001 – Friedmann A, Strietzel FP, Marezki B, Pitaru S, Bernimoulin JP. (2001) Observations on a new collagen barrier membrane in 16 consecutively treated patients. *Clinical and histological observations.* *J Periodontol.* 72(11):1616-1623.
- c. Friedmann *et al.* 2002 – Friedmann A, Strietzel FP, Marezki B, Pitaru S, Bernimoulin JP. (2002) Histological assessment of augmented jaw bone utilizing a new collagen barrier membrane compared to a standard barrier membrane to protect a granular bone substitute Material. *Clin. Oral Impl. Res.* 13, 587–594.
- d. Klinger *et al.* 2010 – Klinger A, Asad R, Shapira L, Zubery Y. (2010) In vivo degradation of collagen barrier membranes exposed to the oral cavity. *Clin. Oral Impl. Res.* 21, 873–876.
- e. Hoang and Mealey 2012 – Hoang TN, Mealey BL. (2012) Histologic comparison of healing after ridge preservation using human demineralized bone matrix putty with one versus two different – sized bone particles. *J Periodontol.* 83:174-181.
- f. Beitlitum *et al.* 2010 – Beitlitum I, Artzi Z, Nemcovsky CE. (2010) Clinical evaluation of particulate allogeneic with and without autogenous bone grafts and resorbable collagen membranes for bone augmentation of atrophic alveolar ridges. *Clin. Oral Impl. Res.* 21: 1242–1250.
- g. Scheyer and McGuire 2014 – Scheyer T. and McGuire MK. (2014) Evaluation of Premature Membrane Exposure and Early Healing in Guided Bone Regeneration of Peri-Implant Dehiscence and Fenestration Defects With a Slowly Resorbing Porcine Collagen Ribose Cross-Linked Membrane: A Consecutive Case Series. *Clinical Advances in Periodontics.* 5: 165-170.

- h. Friedmann et al. 2008 – Friedmann A, Dehnhardt J, Kleber BM, and Bernimoulin, JP. (2008) Cytobiocompatibility of collagen and ePTFE membranes on osteoblast-like cells in vitro. J Biomed Materials Res A 86:935-941.

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Datum Dental Ltd.
MedNet EC-REP III GmbH